

# Non-Acidification Technique for Extracted Chlorophyll *a* Analysis

## Accuracy

Use of the Non-Acidification Optical Kit (P/N 10-040R) will result in more accurate chlorophyll *a* concentration data in all environments. The optical kit is specifically designed to reduce errors due to the background fluorescence from interfering compounds such as chlorophylls *b* and *c*, pheophytin and dissolved organic matter. The improved accuracy is achieved by the following:

**Lamps** – The Blue Mercury Vapor Lamp has narrow lines of emission which optimally excite chlorophyll *a* at 436nm and minimally excite other interfering compounds (see Figure 1).




Figure 1

**Filters** – Narrow band interference filters are selected to only allow the specific excitation (436nm) and emission (680nm) wavelengths of chlorophyll *a* to pass (see Figure 2).




Figure 2

## Approved by EPA

The Non-Acidification Kit is now EPA approved. For Method 445.0 rev. 1.2 or go directly to the EPA webpage at

<http://www.epa.gov/nerlcwww/marinmet.htm>

**NOTE:** the Non-acidification technique is referred to as the “modified fluorometric” technique in Method 445.0.

## Easier

The fluorometric procedure and concentration calculations are simplified due to the absence of the pheophytin correction.

**NOTE:** Pheophytin concentration cannot be determined with the Non-Acidification technique.

## Chlorophyll *b* Interference

The Non-Acidification Kit is not sensitive to high chlorophyll *b* concentrations and is well suited for freshwater and open ocean environments.

## Recent Publications

Tittel, J., B. Zippel, W. Geller. 1998. Relationships between plankton community structure and plankton size distribution in lakes of Northern Germany. *Limnol. Oceanogr.* **43(6)**: 1119-1132.

Welschmeyer, N.A. 1994. Fluorometric analysis of chlorophyll *a* in the presence of chlorophyll *b* and pheopigments. *Limnol. Oceanogr.* **39(8)**: 1985-1992.

Yahel, G., A.F. Post, K. Fabricius, D. Marie, D. Vault, A. Genin. 1998. Phytoplankton distribution and grazing near coral reefs. *Limnol. Oceanogr.* **43(4)**: 551-563.