

AquaFlash™

Handheld Active Fluorometer



Assess Photosynthetic Efficiency of Algae

AquaFlash™ Handheld Active Fluorometer

The AquaFlash™ provides quick and accurate estimations for total chlorophyll and photosynthetic efficiency "health" of algae using *in vivo* fluorescence detection. The sample is saturated with high intensity light to quickly attain the maximum fluorescence state of algae. The fixed saturation pulse is optimized for this package enabling accurate results in less than 15 seconds. Simply insert a 10x10mm square glass or quartz cuvette with your sample and press READ. Results are displayed and automatically logged to be viewed or downloaded at a later time. The AquaFlash has a total data storage capacity of 1,000 measurements which include raw fluorescence values along with calculated estimates for user reference.

Ordering Information

INSTRUMENT	PART NUMBER
AquaFlash™ Handheld Active Fluorometer	8600-000
ACCESSORIES	PART NUMBER
Wrist Strap	030-8500
AquaFlash™ Calibration Solution	8600-225

Highlights of the AquaFlash™

- Fast
 - Results in <15 seconds
 - Total Chlorophyll (µg/L)
 - Photosynthetic Efficiency (yield)
 - Quick Calibration Check
- Portable
 - Handheld, Battery Powered
 - Highly Durable Case
 - Internal Data Logging

AquaFlash™ Specifications

MDL	0.3 µg/L
Linear Range	0.3-100 µg/L
Linearity	0.99R ²
Weight in Air	13.9 oz; 0.4 kg
Size	1.75" x 3.5" x 7.25"; 4.45 cm x 8.9 cm x 18.4 cm
Warm-up Time	5 seconds
Case	IP 67 standard; dustproof/waterproof
Temperature	41-104°F; 5-40°C
Power	4 AAA batteries (standard or rechargeable)
Max Data Capacity	1,000 measurements
Data Output	ASCII

AquaFlash limits were determined using *P. micans* and *Tetraselmis sp.*; limits may vary for other algal species

Contact Us

Toll-Free : 1.877.316.8049

Email : sales@turnerdesigns.com

Address:

Phone : 408.749.0994

Web : www.turnerdesigns.com

1995 N. 1st Street

Fax : 408.749.0998

San Jose, CA 95112