

**Trilogy Power Supply Kit (P/N 7200-941)**

This kit is included as an accessory with the Trilogy. It contains interchangeable clips for US, UK, European and Australian wall sockets and voltage stabilization for 10-12V power sources.



**Trilogy Data Cable (P/N 021-0700)**

This kit is included as an accessory with the Trilogy. It is a DB9 male to female cable used for data transfer from the Trilogy to a PC.



**If your computer only has USB ports we recommend purchasing a USB to serial adapter cable for use with the Trilogy Data Cable. <http://www.belkin.com/> search for Part Number: F5U409v1.**

